

Jabra Elite 7 Active - Mint

How do I achieve the best result with MyFit?

When you start MyFit, a multitone is played in the earbuds to detect if the seal of each earbud can be improved. Therefore, the most important thing is not to touch your earbuds when MyFit is in progress.

If MyFit detects that the seal in your earbuds can be improved, you will be asked to adjust the position of the earbuds, change the EarGel size, or do both. Make sure you run MyFit with all three sizes of EarGels until you achieve a good fit.

Note

If MyFit cannot identify the best seal for you, we recommend that you use the EarGels that you feel give the best combination of comfort and sound quality.

Applicable products



Jabra Evolve2 Buds -
USB-C MS - Wireless



Jabra Evolve2 Buds -
USB-A UC - Wireless



Jabra Evolve2 Buds -
USB-C UC



Jabra Evolve2 Buds -
USB-A MS - Wireless



Jabra Evolve2 Buds -
USB-C MS



Jabra Evolve2 Buds -
USB-A UC



Jabra Evolve2 Buds -
USB-A MS



Jabra Evolve2 Buds -
USB-C UC - Wireless



Jabra Elite 85t -
Titanium Black



Jabra Elite 85t - Grey



Jabra Elite 85t - Grey
with 2 wireless



Jabra Elite 7 Sport -
Black



Jabra Elite 7 Pro -
Black



Jabra Elite 7 Pro -
Gold Beige



Jabra Elite 7 Active -
Navy



Jabra Elite 7 Active -
Black



Jabra Elite 85t - Black



**Jabra Elite 85t -
Copper Black**



**Jabra Elite 7 Pro -
Titanium Black**



**Jabra Elite 85t - Gold
Beige**



**Jabra Elite 7 Active -
Mint**